

BPET Series - Extraction of damaged probes



In some cases, probes can break inside of a receptacle. It can be challenging to remove damaged barrels especially in case of threaded probes without the correct tool. We can eliminate this problem with our BPET series. BITS have special left-handed threaded tips to make extraction simple. Just insert the tool and press down a little while turning it counter-clockwise, then remove the damaged barrel with ease.

BPET Series:

- BPET-102 - Extraction of damaged probes with $\varnothing 1,02\text{mm}$
- BPET-135 - Extraction of damaged probes with $\varnothing 1,37\text{mm}$
- BPET-165 - Extraction of damaged probes with $\varnothing 1,66\text{mm}$
- EX-14 - Extraction of damaged probes with $\varnothing 2,65\text{mm}$

Replacement BIT (BIT-BPET-xxx) and handle (H-2.6) can be ordered separately.